



*Everything For Your Implant
Practice But The Implants®*

2024 REGENERATIVE CATALOG





Dr. John Russo
Periodontist
Sarasota, FL

"Salvin Dental is a great resource for my practice and especially for the doctors who take my courses. Salvin has what you need whether you're getting started or doing advanced procedures."



Dr. Michael Pruett
General Dentist
Augusta, GA

"I love the way that Salvin does business. Their willingness to go the extra mile and their commitment to serving makes it a pleasure to deal with them."



Dr. Nick Shumaker
Periodontist
Fort Collins, CO

"When it comes to quality products and great customer service, nobody does it better than Salvin!"



Dr. Fred Sakamoto
Periodontist
Blacklick, OH

"Salvin Dental runs a first class operation with excellent customer service and products. We consider Salvin as part of our team with periodontal, implant and related surgical cases."




Dr. Kanyon Keeney
Oral Surgeon
Mechanicsville, VA

"Salvin Dental has one of the most professional sales teams in the dental industry. They know that my time is valuable and always deliver more than they promise."



Dr. Lauralee Nygaard
Periodontist
Spokane Valley, WA

"When Salvin says everything for the implant practice but the implants, they mean it! Salvin is our prime resource for expertise and incredible service. They are a doctor's true friend in our industry."



Dr. Etienne Lacrampe
Periodontist
Walnut Creek, CAI

"I use Salvin for all my biologics. There is no one I trust more to provide the best results for both my patients and my referring dentists."



Dr. Richard Martin
Oral Surgeon
Lewisville, TX

"When you work with Salvin Dental, you do not feel like you are working with a company, you feel like you are working with your family. Their attention to detail and reliability is unmatched in the dental field."



Dr. Eric Hamrick
Periodontist
Greenville, SC

"For over twenty-five years, Salvin has been our practice's go-to company for surgical supplies and instruments. Their service, vast inventory, and prices are comparable to none."



Dr. Avi Schetritt
Periodontist
Miami, FL

"Great products, great service, and great people! Their professionalism is what sets Salvin apart. They deliver every time!"



Dr. Bradley Ross
Periodontist
Miami, FL

"Salvin provides top quality surgical armamentarium and biomaterials, but their people truly make this company an essential asset to my practice. Their knowledge, integrity, and customer service are second to none."



Dr. Matthew Fien
Periodontist
Deerfield Beach, FL

"I am always impressed with the professional service and knowledge from the team at Salvin Dental. They make it easy for me to recommend them to my colleagues."



Dr. Steve Wallace
Periodontist
Wilmington, NC

"I turn to Salvin when I need answers to my questions. Their knowledge and helpfulness make them a pleasure to do business with."



Dr. Justin Zalewsky
Periodontist
Alexandria, VA

"Salvin has been a business partner for years, because they truly do partner with us. Their products and services are unmatched. Their people care about our practice as much as we do."



Dr. Nancy Newhouse
Periodontist
Independence, MO

"Nobody does it better! Salvin is just what they say they are: 'Everything For Your Implant Practice But The Implants' all wrapped up in the best service in the industry."



Dr. David Lee Hill
Oral Surgeon
Chapel Hill, NC

"I can't think of a single source that offers a more comprehensive quality line of products to support my implant practice than Salvin Dental. They always have what I need."



Dr. James Woodyard
Periodontist
Newburgh, IN

"Knowing that I can expect the right products and great service is why I choose to do business with Salvin. Dealing with them is a first-class experience every time I call."



Dr. James Rutkowski
General Dentist
Clarion, PA

"Honesty and integrity are values that come to mind when working with Salvin Dental. I know that I'm going to be taken care of every time I call."



Dr. Israel Puterman
Periodontist
Chevy Chase, MD

"Getting the Salvin catalog is like being a kid again, flipping through the toy catalog and seeing all the things that I want. The quality of the products is only surpassed by the quality of the people. They truly want your success above all else."



Dr. Paul Fugazzotto
Periodontist
Milton, MA

"Salvin has, in my experience with them over decades, consistently demonstrated a degree of honesty and openness, and a genuine desire to help the clinician, that is unsurpassed in the dental community."



Dr. Jason Souyias
Periodontist
Port Huron, MI

"The team at Salvin is amazing! They know customer service better than anyone in dental. That is why they are so easy to work with."

SALVIN[®] Dental Specialties

SALVIN[®] Regenerative

Everything For Your Implant Practice But The Implants[®]

2024 REGENERATIVE PRODUCTS CATALOG

*We're in business to make your
lives easier and improve your
surgical outcomes.*

Call us, we'll take great care of you!

—Bob Salvin



Please See Our 2024 Instrument Catalog For:

- *Instruments*
- *Vital Sign Monitors*
- *Tacks & Screws*
- *Sutures*
- *Motors & Handpieces*
- *Infection Control*
- *O/R Furniture*
- *Torque Wrenches & More*

800-535-6566 [Toll Free U.S. & Canada • www.salvin.com



Salvin Dental Specialties, Inc.
3450 Latrobe Drive
Charlotte, NC 28211 USA
Ph: 704-442-5400 Fax: 704-442-5424
Email: orders@salvin.com
Hours: M-Th: 8:30-7:00 F: 8:30-6:00
(U.S. Eastern Standard Time)

Quality Management System
ISO 13485:2016 Certified

*See It Demonstrated!
Look For This Icon Throughout
Our Catalog*



Special Thanks To Susan Salvin,
Our Cover Model.

Portrait Photo By Bob Salvin.

Stylized Photo Illustration By
Karen Colella, Senior Art Director.

Trademarks: "Salvin[®]", "Salvin Dental Specialties[®] Logo", "Everything For Your Implant Practice But The Implants[®]", "Salvin Regenerative[®]", "Salvin Regenerative[®] Logo", AlloSculpt[®], Salvin AmnioExcite[®], CytoSurg[®], OraCOTE[®], OraFoam[®], OraMEM[®], OraMEM-Sustained[®], OraPLUG[®], OraTAPE[®], Orastruct[®], OsseoFiber[®], OsseoFiber-Plus[®], Renovix-Plus[®], SalvinOss[®], SalvinOss[®] Collagen, and Xymphony[®] are registered U.S. and international trademarks of Salvin Dental Specialties, Inc.

Salvin Amplifi[™] is trademarks of Salvin Dental Specialties, Inc. Other product brand names may be trademarks or registered trademarks of their respective owners and are mentioned for reference purposes only.

CONTENTS

5-6

AmnioExcite Full Thickness Placental Membrane

Trophoblast Layer Contains Significant Biologically Relevant Proteins

4X Thicker Than Amnion Or Chorion For Much Better Handling



7

OsseFiber Demin Cortical Fibers

Mix With Mineralized Graft Materials For GBR

Handles Like Putty - 100% Bone - No Carrier



8

OsseFiber-PLUS Demin Cortical Fibers + Mineralized Cancellous Chips

For GBR, Ridge Augmentation, Sinus Lifts & Bony Defects

Handles Like Putty - 100% Bone - No Carrier



9-12

LifeNet Health Allograft Bone

A Complete Selection Of LifeNet Health Allografts



13-14

SalvinOss Xenograft Particulate

For Sinus Lifts, Ridge Augmentation & GBR

Maintains Space During Bone Remodeling

High Porosity




15-16

SalvinOss Collagen Xenograft + Collagen Plug

For Socket Preservation & Perio Defects

Retains Shape After Hydrating With Blood Or Saline



17

Salvin Amplifi Long Term Membrane

Resorbs in 6-9 Months For Long-Term Barrier Functionability

Significant Tensile Strength



18

Renovix-Plus Resorbable Extra Cellular Matrix

Maintains Barrier Protection For Up To 6 Months

Our Most Versatile Resorbable Membrane



19

OraMEM Cross Linked Bovine Collagen

Tuck & Place In Buccal Defects

Retains Shape After Hydrating

Resorbs In 4-8 Weeks



19-20

OraMEM-Sustained Cross Linked Bovine Collagen

For Sinus Lifts & Larger GBR Procedures

Retains Shape After Hydrating

Resorbs In 18 Weeks



21

Xymphony Drapeable Porcine Pericardium

33 Week Barrier Function

Significant Tensile Strength



22

CytoSurg Non-Resorbable PTFE Membrane

For Socket Preservation When
Primary Closure Cannot Be Achieved
Remove In 2-3 Weeks



23-24

OraCELL Decellularized Dermis Allograft

For GTR & Soft Tissue Grafts
Prehydrated - Ready To Use Right Out
Of The Package
Store At Room Temperature
No Refrigeration Required



25

OraPLUG Absorbable Collagen Plug

To Stabilize Blood Clots & Preserve
Graft Material In Socket Preservation



25

OraTAPE Absorbable Collagen Wound Dressing

To Repair Sinus Membrane Perfs &
Stabilize Blood Clots In Oral Wounds



26

FULFILL 80% Calcium Phosphate Granules Combined With 20% Bovine Collagen

Mimics Natural Bone In Resorption &
Remodeling Profile



26

OraCOTE Absorbable Collagen Wound Dressing

For Palatal Donor Sites & Stabilize
Blood Clots In Oral Wounds
3X Thicker Than OraTAPE®



27

AlloSculpt-3D Demin Cortical With Min Cortical + Cancellous Chips

For GBR, Ridge Augmentation,
Sinus Lifts & Bony Defects
Handles Like PlayDoh™ When
Hydrated

Hydrogel Carrier Helps
Keep Graft Material In Place



28

Accell DBM Putty In A Syringe

For GBR, Ridge Augmentation,
Sinus Lifts & Bony Defects

Each Lot Tested For
Osteoinductive Potential



29

DynaBlast Demin Cortical + Min Cancellous Chips Paste or Putty

For GBR, Ridge Augmentation,
Sinus Lifts & Bony Defects

Carrier Thickens At Body Temp To
Keep Graft In Place



29

DynaGraft Demin Cortical Paste or Putty

For Osseous Defects
Carrier Thickens At Body Temp
To Keep Graft In Place



30

Orastruct 80/20 TCP/HA Graft Material

Alternative To Allografts Or Xenografts
For GBR

Slow Resorbing Osteoconductive
Space Maintaining



30

OraFOAM Hemostatic Absorbable Collagen Sponge

Use In Extraction Sockets


Controls Bleeding In 2 Minutes



31-34

PTFE Sutures, Sutures & Wound Closing






Full-Thickness Placental Membrane

- 4X Thicker Than An Amnion Or Chorion Membrane Making It Much Easier To Handle
- Does Not Fold In Upon Itself When Rehydrated

Structure of Salvin AmnioExcite® Full-Thickness Placental Membrane



*AmnioExcite® Includes Both Amnion And Chorion PLUS The Trophoblast Layer
The Trophoblast Layer Contains More Than 50% Of
The Biologically Relevant Proteins In The Placental Membrane (1)*


SALVIN AMNIO EXCITE Amnion / Chorion Membrane

Size	Product #	Fee	7 or More
10mm x 10mm	#AMNIO-MEM-10X10	97.00	86.00
15mm x 20mm	#AMNIO-MEM-15X20	184.00	170.00
20mm x 30mm	#AMNIO-MEM-20X30	283.00	261.00



1. LifeNet Health®. TR-004-2020 Characterization of the Amnion, Chorion, and Trophoblast Layers of Decellularized and Freeze-Dried Placental Membrane. 2020.





AmnioExcite® Case Studies

Case Study #1 - Socket Preservation Using AmnioExcite® Full-Thickness Placental Membrane


Case & Photos: Dr. Nick Shumaker



Surgical Extraction



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*



*3 Week Post-Op –
Excellent Tissue Healing*

Case Study #2- Cemental Tear/Osseous Defect Repair Using AmnioExcite®

Case & Photos: Dr. Nick Shumaker



*Cemental Tear /
Osseous Defect*



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*




*6 Month - No Recession
Excellent Aesthetic Result*



“As a user of amnion membranes for over 15 years, I have found AmnioExcite® from Salvin to be a great addition to my practice. AmnioExcite® is much easier to handle and place, but still adapts to adjacent structures such as bone, tooth roots or dental implants. It has the perfect balance of adaptation and handling that has been missing in amnion membranes, and I’m very pleased with the results.”

—**Dr. Nick Shumaker**
Periodontist, Fort Collins, CO



Demineralized Cortical Fiber

- *Demineralized Cortical Bone Has The Potential To Stimulate New Bone Growth*
- *100% Bone - No Carrier*
- *Handles Like Putty When Hydrated*
- *Mixes Well With Mineralized Graft Materials If Needed*

100% Bone – No Carrier – Handles Like Putty When Hydrated Facilitates Cellular Proliferation & Cellular Attachment ⁽¹⁾



SALVIN OsseoFiber™ Demineralized Cortical Fibers



Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-0.5CC	153.00	136.70
1.0 cc	#OFIBER-1.0CC	261.00	233.80
2.5 cc	#OFIBER-2.5CC	512.00	464.60


No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 3 Years



**Donor Tissue From LifeNet Health® Tissue Processor
Each Lot Is Tested In Vivo For Osteoinductive Potential**

1. Rodriguez RU, Kemper N, Breathwaite E, et al. Demineralized bone matrix fibers formable as general and custom 3D printed mold-based implants for promoting bone regeneration. *Biofabrication*. 2016;8(3):035007





Demineralized Cortical Fiber PLUS Mineralized Cancellous Chips

- *Demineralized Cortical Bone Has The Potential To Stimulate New Bone Growth*
- *Mineralized Cancellous Chips Provide Scaffold For Bone Cell Migration*
- *Use For GBR Procedures, Ridge Augmentation, Sinus Lifts & Bony Defects*



100% Bone – No Carrier – Handles Like Putty When Hydrated Facilitates Cellular Proliferation & Cellular Attachment ⁽¹⁾

SALVIN OsseoFiber-PLUS™ Demineralized Cortical Fibers PLUS Mineralized Cancellous Chips



Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-PLUS-0.5CC	159.00	141.70
1.0 cc	#OFIBER-PLUS-1.0CC	273.00	245.80
2.5 cc	#OFIBER-PLUS-2.5CC	514.00	465.50


No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 3 Years



Moldable & Easily Shaped To Pack Into Any Defect

Case Study - Guided Bone Regeneration Using OsseoFiber-Plus™


Case & Photos: Dr. Paul Rosen




Post Surgical Extraction – Osseous Defect



OsseoFiber-Plus™ Has Putty-Like Handling



Successful GBR 6 Month Post-Op



Implant Placement Prior To Restoration



“The innovation of making DFDBA in a fiber form has been a tremendous enhancement to this graft. The increase in surface area makes this osseoinductive clot stabilizer a must for one's regenerative armamentarium. Importantly, this fiber form of DFDBA transforms it into being much more user friendly.”


–Dr. Paul Rosen

Periodontist, Yardley, PA & New York, NY


**Donor Tissue From LifeNet Health® Tissue Processor
Each Lot Is Tested In Vivo For Osteoinductive Potential**

1. Rodriguez RU, Kemper N, Breathwaite E, et al. Demineralized bone matrix fibers formable as general and custom 3D printed mold-based implants for promoting bone regeneration. *Biofabrication*. 2016;8(3):035007





Over 10 Million Bio-Implants Have Been Processed Using LifeNet Health's® Allowash Technology, With No Disease Transmission



LifeNet Health® Is The World's Largest Full Service Tissue Processor / Organ Procurement Organization


Of Total Donors Screened For Transplantation, Only 3% Are Accepted By LifeNet Health®

All Of LifeNet Health's® Donor Tissue Is Recovered In The United States

Specific Screening Criteria*	LifeNet Health® Saving Lives. Restoring Health. Giving Hope.	AATB (American Association of Tissue Banks)	FDA (U.S. Food & Drug Administration)
Encephalitis (clinically active)	●	●	●
Gonorrhoea (within the past year)	●	●	●
Hepatitis B	●	●	●
Hepatitis C	●	●	●
HIV 1/2	●	●	●
Malaria	●	●	●
Meningitis (clinically active)	●	●	●
Sepsis	●	●	●
Syphilis	●	●	●
Transmission spongiform encephalopathy (TSE)	●	●	●
Vaccinia	●	●	●
West Nile Virus (WNV, 4 month exclusion from the time of onset)	●	●	●
Guillain-Barre syndrome (clinically active)	●	●	●
Multiple sclerosis	●		●
Leprosy (Hansen's disease)	●	●	
Clinically significant metabolic bone disease	●	●	
Polyarteritis nodosa	●	●	
Rabies	●	●	
Rheumatoid arthritis (musculoskeletal deferral)	●	●	
Sarcoidosis	●	●	
Systemic lupus erythematosus	●	●	
Systemic mycosis	●	●	
Tuberculosis (clinically active or treated in the past 5 years)	●	●	
Ankylosing spondylitis	●	●	
Autoimmune vasculitis if type is not known	●	●	
Cancer	●	●	
Chagas disease	●	●	
Epstein Barr virus (clinically active)	●	●	
Clostridium difficile infection (defer if occurred during present hospitalization)	●	●	
Endocarditis (clinically active)	●	●	
Herpes (any evidence of active herpes infection is a deferral)	●	●	
Illicit drug use (see Illegal drugs)	●	●	
Peritonitis (clinically active)	●	●	
Rheumatic fever (defer all tissues if clinically active, defer heart valves for any history)	●	●	
Poliomyelitis (clinically active or untreated)	●	●	
Mixed connective tissue disease	●		
Cold agglutinin disease	●		
Pyelonephritis (clinically active)	●		
Reactive arthritis (must be asymptomatic for 6 months for musculoskeletal deferral)	●		
Wegener's granulomatosis	●		

*These Screening Criteria Are Subject To Change

LifeNet Health®, The LifeNet Health® Logo, OrACELL®, OraGRAFT® & Preservon® Are Registered Trademarks Of LifeNet Health®.



LifeNet Health® OraGRAFT®

Mineralized 50/50 Cancellous/Cortical Mix **MOST POPULAR**



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#CAN-CORT-.25	65.00	58.00
	0.50 cc	#CAN-CORT-.50	96.00	85.60
	1.00 cc	#CAN-CORT-1.00	129.00	115.80
	2.00 cc	#CAN-CORT-2.00	186.00	167.40

Mineralized Cancellous **MOST POPULAR**




Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CAN-.25	86.00	76.80
	0.50 cc	#MIN-CAN-.50	99.00	87.70
	1.00 cc	#MIN-CAN-1.00	134.00	119.40
	2.00 cc	#MIN-CAN-2.00	206.00	184.60
1000-2000 microns	0.50 cc	#MIN-CAN-L-.50	99.00	87.70
	1.00 cc	#MIN-CAN-L-1.00	134.00	119.40
	2.00 cc	#MIN-CAN-L-2.00	206.00	184.60

Mineralized Cortical



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CORT-.25	47.00	40.80
	0.50 cc	#MIN-CORT-.50	65.00	57.40
	0.70 cc	#MIN-CORT-.75	81.00	71.70
	1.20 cc	#MIN-CORT-1.25	95.00	84.50
	2.50 cc	#MIN-CORT-2.50	140.00	124.80
250-710 microns	0.50 cc	#MIN-CORT-SM-.50	65.00	57.40
	0.70 cc	#MIN-CORT-SM-.75	81.00	71.70
	1.20 cc	#MIN-CORT-SM-1.25	95.00	84.50
	2.50 cc	#MIN-CORT-SM-2.50	140.00	124.80





LifeNet Health® OraGRAFT®

70/30 Mineralized/Demineralized Cortical Mix **MOST POPULAR**



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.50 cc	#MD-7030-.50	92.00	82.60
	1.00 cc	#MD-7030-1.00	134.00	120.50
	2.00 cc	#MD-7030-2.00	189.00	170.70

Demineralized Cortical



Particle Size	Size	Product #	Fee
250-1000 microns	0.25 cc	#DEMIN-CORT-.25	53.00
	0.50 cc	#DEMIN-CORT-.50	71.00
	0.70 cc	#DEMIN-CORT-.75	90.00
	1.20 cc	#DEMIN-CORT-1.25	105.00
	2.50 cc	#DEMIN-CORT-2.50	155.00

Demineralized Cancellous




Particle Size	Size	Product #	Fee
1-8 mm Crushed	5.00 cc	#DEMIN-CAN-CR-5	324.00

Graft Site Suggested Volume Guide

Periodontal Defects		Extraction Sites		Sinus Grafts	
Central Incisor	0.5 cc	Central Incisor	0.7-1.0 cc	Each Sinus	2.5-5.0 cc
Lateral Incisor	0.5 cc	Lateral Incisor	0.5 cc		
Canine	0.5 cc	Canine	0.7-1.0 cc		
1st & 2nd Premolar	0.5 cc	1st & 2nd Premolar	0.7-1.0 cc		
1st & 2nd Molar	0.7-1.0 cc	1st & 2nd Molar	1.2-2.0 cc		
3rd Molar	0.7-1.0 cc	3rd Molar	1.2-2.0 cc		

For optimal safety, all LifeNet Health® allografts distributed by Salvin Dental are gamma irradiated to achieve a Sterility Assurance Level (SAL) 10⁻⁶ without compromising the inherent biomechanical or biochemical properties of the graft. Typical Shelf Life: 4 Years. Fees are subject to change. Please check current fees at time of placing order. LifeNet Health® Canadian Registration: CTO #100038



LifeNet Health® OraGRAFT®

Pre-Hydrated Cortical Plate

Ready To Use Out Of The Package

- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature
- Preservon-Treated Grafts Have Been Shown To Have Greater Flexibility Than Freeze-Dried Grafts ⁽¹⁾
- Ideal For “Shell” Technique – Secure Cortical Plate To Host Bone With Fixation Screws To Create A Biologic Container



Plate	Size	Product #	Fee
Cortical Plate	30 mm x 15 mm x 1 mm	#CORT-PLATE-30X15X1	510.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 5 Years

Pre-Hydrated Ilium Block

Ready To Use Out Of The Package

- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Ilium Block - Small	15 mm x 20 mm	#ILIUM-PRES-SMALL	571.00
Ilium Block - Medium	15 mm x 30 mm	#ILIUM-PRES-MEDIUM	648.00
Ilium Block - Large	20 mm x 40 mm	#ILIUM-PRES-LARGE	864.00

Pre-Hydrated Unicortical Block MOST POPULAR

Ready To Use Out Of The Package


- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Unicortical Block - Small	10 mm x 15 mm x 8 mm	#ILIUM-UNICORTICAL-SMALL	329.00
Unicortical Block - Large	10 mm x 25 mm x 8 mm	#ILIUM-UNICORTICAL-LARGE	399.00

Case Study - GBR Using LifeNet Health® Unicortical Block & Renovix-PLUS® Membrane


Case & Photos: Dr. Shaun Young




Large Osseous Defect



LifeNet Health®
Unicortical Block



Renovix-Plus® –
Barrier Membrane




Implant Placed In
Regenerated Bone



“Salvin continues to raise the bar with their regenerative products. The allogeneic blocks have elevated my abilities and given me additional options that I can depend on for ridge augmentation. The Renovix-PLUS® collagen membrane is one of the best on the market due to its versatility as a stiffer membrane, but also very moldable – it’s really the best of both worlds.”


–Dr. Shaun Young
Oral Surgeon, Tampa, FL

1. LifeNet Health. TR 0216 Three-Point Bend Testing of Glycerol-Treated, Freeze-Dried, and Frozen Bone



SalvinOss® Mineralized Cancellous Xenograft

- Indicated For Sinus Lifts, Ridge Augmentation & Guided Bone Regeneration
- Maintains Space & Facilitates Osteoconduction
- 0.25 - 1.0mm High Porosity Particles Become An Integral Part Of Newly Formed Bone
- Radiopaque



Exceptional Handling Characteristics


SALVIN OSS Xenograft Bone Graft Material



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.50 CC (0.25G)	#SALVINOSS-0.50	96.00	86.60	75.50
	1.00 CC (0.5G)	#SALVINOSS-1.00	144.00	130.40	117.70
	2.50 CC (1.25G)	#SALVINOSS-2.50	263.00	239.50	216.60

Download Our White Paper At: www.salvin.com/Downloads/Guides-Studies.html






SalvinOss® Case Studies

Case Study - Socket Preservation With SalvinOss®


Case & Photos: Dr. Michele Holzinger



Pre-Op



Post-Extraction –
4 Months Healing



Post-Implant –
3 Months Healing




"I trust Salvin Regenerative because they help me produce results that are consistent and reliable for my patients. When I re-enter a site grafted with SalvinOss®, I know exactly what to expect because the results are **predictable**."


–Dr. Michele Holzinger
Periodontist, Middletown, CT

Case Study - Ridge Augmentation Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Mitchel Godat



Severely Pneumatized Sinus Cavity



12 Month Post-Op Radiograph




"I have been using SalvinOss® for many years to increase bone volume for space preservation prior to implant placement. SalvinOss® is excellent for sinus lifts, ridge augmentation, and other GBR / GTR procedures with very predictable long-term results."


–Dr. Mitchel Godat
Periodontist, Memphis, TN

Case Study - Guided Tissue Regeneration Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Grant King



9-10 mm
Periodontal Defect



8 Month Post-Op –
Significant Pocket Fill



"SalvinOss® handles extremely well, especially when I mix it with blood, making it easy for me use for GTR and GBR. It is a slower resorbing material, holding space and volume more than using allograft alone. Most importantly, SalvinOss® gives me predictable results every time."

–Dr. Grant King
Periodontist, Memphis, TN

SalvinOss® Collagen On Next Page →

SalvinOss® Collagen Xenograft Bone + Bovine Collagen

- 90% SalvinOss® Xenograft Bone Graft Material PLUS 10% Bovine Collagen
- For Socket Preservation, Sinus Lifts, And Periodontal / Furcation Defects
- Plug Shape Helps Make Socket Preservation Fast & Simple
- Retains Shape After Hydrating



7mm

10mm

100mg

Ideal For:

- Anterior Socket Grafting
- Periodontal Defects
- Osseous Defects Caused By Peri-Implantitis

10mm


10mm

250mg


Ideal For:

- Posterior Socket Grafting
- Ridge Augmentation
- Sinus Lifts

Plugs Shown Actual Size



SALVIN(OSS) Collagen Xenograft Bone PLUS Bovine Collagen



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	7mm (Diam) x 10mm	#SALVINOSS-COLL-100MG	113.00	106.50	100.50
	10mm (Diam) x 10mm	#SALVINOSS-COLL-250MG	200.00	192.00	182.00

Download Our White Paper At:
www.salvin.com/Downloads/Guides-Studies.html

Scan For


60th Second Product Video™

Scan For

60th Second Product Video™

Scan For

SalvinOss® Collagen Technique Guide



SalvinOss® Collagen Case Studies

Case Study #1 - Socket Preservation Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg



Plug Shape For Easy Socket Preservation And Crestal Sinus Lift

Case Study #2 - Furcation Defect Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg




SalvinOss® Collagen Helps When Graft Containment Is Problematic



"SalvinOss® Collagen has become a staple of my practice for regeneration as well as augmentation. The ease of handling coupled with the volumetric stability produces predictability in grafting situations around both teeth and implants alike."

-Dr. Tyler Borg
Periodontist, Northglenn, CO



Membranes Shown Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm

Resorbs In 6-9 Months For Long-Term Barrier Functionability

- Collagen-Fiber Matrix With Intact Intrinsic Crosslinks For Mechanical Strength
- <0.01 Micron Pore Size Allows Nutrients & Growth Factors While Assuring Cell Barrier Function
- Chemically Engineered To Control Stability & Minimize Inflammatory Response
- Significant Tensile Strength For Ease Of Tacking & Suturing

SALVIN Amplifi Long Term Membrane		Resorbs In 6 To 9 Months		
Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#AMPLIFI-15X20	168.00	149.70	138.60
20mm x 30mm	#AMPLIFI-20X30	210.00	187.40	174.20
30mm x 40mm	#AMPLIFI-30X40	305.00	278.80	259.70




Case Study - Amplifi™ Long-Term Membrane


Case & Photos: Dr. Bradley Ross



Surgical Extraction



Bone Graft Material In Defect



Amplifi™ Sutured In Place




New Bone Formation – Ready For Implant



“The Amplifi™ membrane from Salvin has become a go to biomaterial for me for ridge augmentation procedures. It is extremely easy to manipulate and use, is very resilient, and most importantly it has a long resorption time. The results I have seen with the Amplifi™ membrane have been spectacular.”

–Dr. Bradley Ross
Periodontist, Miami, FL

*Data On File



Our Most Versatile Resorbable Barrier Membrane

- For Bone Regeneration & Soft Tissue Esthetics
- Non-Cross Linked Extracellular Matrix Perfectly Conforms To The Defect When Hydrated
- Contains Type I, III, & IV Collagen
- Excellent Tensile Strength & Tear Resistance
- Trim & Place Either Wet Or Dry
- Resorbs Within 6 Months

Case Study - Renovix-Plus® Guided Healing Collagen

Case & Photos: Dr. Lauralee Nygaard




Aesthetic Zone
Surgical Extraction



Renovix-Plus®
Over Graft Material



Soft Tissue Aesthetics
Preserved



Implant Placed



"Renovix-Plus® has become the membrane I use most for predictable ridge preservation. I use it for ease of handling and preservation of soft tissue architecture. And the bonus of working with team Salvin!"

–Dr. Lauralee Nygaard
Periodontist, Spokane Valley, WA

SALVIN® Renovix-Plus® Guided Healing Collagen Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#RENOVIX-PLUS-15X25	177.00	159.40	148.60
20mm x 30mm	#RENOVIX-PLUS-20X30	231.00	210.50	196.40
30mm x 40mm	#RENOVIX-PLUS-30X40	301.00	268.70	255.50





Resorbable Collagen Membrane

- For Use In Buccal Defects & Smaller GBR Procedures
- Cross Linked Type 1 Bovine Collagen
- Retains Shape & Memory For Easy Placement
- Resorbs In 4-8 Weeks



Resorbs In 4 To 8 Weeks

Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#ORAMEM-15X20	169.00	151.60	140.40
20mm x 30mm	#ORAMEM-20X30	213.00	191.30	178.00
30mm x 40mm	#ORAMEM-30X40	307.00	279.70	260.50

Membranes Shown Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm



Sustained Duration Resorbable Collagen Membrane

- For Use In Sinus Lifts & Larger GBR Procedures
- Cross Linked Type 1 Bovine Collagen
- Retains Shape & Memory For Easy Placement
- Thicker Than Regular OraMEM®
- Resorbs In 18 Weeks




Resorbs In 18 Weeks

Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#ORAMEM-SUS-15X20	196.00	175.40	164.50
20mm x 30mm	#ORAMEM-SUS-20X30	260.00	234.60	219.30
30mm x 40mm	#ORAMEM-SUS-30X40	342.00	308.30	290.40

Scan For

60 Second Product Video™






OraMEM-Sustained® Case Studies

OraMEM-Sustained® Case Study #1 - Socket Preservation / Buccal Wall Repair

Case & Photos: Dr. Matthew Fien




Bone Sounding Revealed Buccal Wall Dehiscence



Socket Grafted With LifeNet Health® 70/30 Bone – OraMEM-Sustained® Tucked Into Buccal Pouch



4 Month Post-Op Beautiful Soft Tissue Aesthetics



4 Month Post-Op – Complete Regeneration of Buccal Wall




“Following atraumatic extraction and debridement of tooth #8, bone sounding revealed a moderate buccal dehiscence. The socket was grafted with particulate bone graft (LifeNet Health® 70/30). A 15x20 OraMEM-Sustained® membrane was trimmed and inserted into a buccal pouch. A resorbable monofilament suture was used to stabilize the site with intentional secondary intention healing. After 4 months, implant surgery was performed and complete regeneration of the buccal wall defect was noted.”

–Dr. Matthew Fien


Periodontist, Fort Lauderdale, FL

OraMEM-Sustained® Case Study #2 - Socket Preservation / Buccal Wall Repair


Case & Photos: Dr. Matthew Fien




Extraction Of Tooth #3 & #4 With Buccal Wall Defect



Grafted With LifeNet Health® 70/30 Bone & OraMEM-Sustained®



4.5 Month Post-Op – Complete Regeneration Of Osseous Defects




Two Implants Placed In New Host Bone



“Following extraction of tooth #3 and #4, a moderate buccal plate dehiscence was noted. The sockets and buccal defects were grafted with particulate bone graft (LifeNet Health® 70/30). An OraMEM-Sustained® membrane was trimmed and adapted over the buccal occlusal line angle. A combination of a monofilament resorbable suture and PTFE was used to stabilize the site without primary closure. After 4.5 months of healing, the osseous defects were completely regenerated and two implants placed.”


–Dr. Matthew Fien

Periodontist, Fort Lauderdale, FL



Our Most Drapeable, "No Shape Memory" Barrier Membrane

- Porcine Pericardium Maintains Barrier Function Up To 33 Weeks
- Significant Tensile Strength - Suture Or Tack With Confidence





SALVIN X YMPHONY[®] Contour Adapting Resorbable Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#XYMPHONY-15X25	160.00	142.40	131.70
20mm x 30mm	#XYMPHONY-20X30	212.00	189.70	176.40
30mm x 40mm	#XYMPHONY-30X40	290.00	260.70	243.80



1. In vivo animal study data on file.
 2. Data on file.





Non-Resorbable PTFE Membrane

- For Socket Preservation When Primary Closure Cannot Be Achieved
- Sufficiently Firm To Protect The Defect – Supple Enough To Drape Over Surrounding Tissue
- Nominal Thickness Of 250 Microns
- Needs To Be Removed In 2-3 Weeks

Space-Making Barrier In The Treatment Of Periodontal Defects

Textured Top Surface Provides Scaffold For Soft Tissue

Smooth Bottom Surface Against Bone For Easy Retrieval



Pore Size Less Than 0.3 Microns



Non-Resorbable PTFE Membrane 1 Membrane Per Box

Size	Product #	Price	10 or More	20 or More
12mm x 24mm	#CYTOSURG-12X24	50.00/Membrane	42.70/Membrane	37.10/Membrane
25mm x 30mm	#CYTOSURG-25X30	86.00/Membrane	75.80/Membrane	66.70/Membrane





Decellularized Dermis Allograft

- Excellent For GTR/GBR & Oral Soft Tissue Correction
- Terminally Sterilized For Increased Patient Safety
- Handles Like Autogenous Tissue




**Prehydrated & Ready To Use
Right Out Of The Package**




**No Refrigeration Required –
Store At Ambient Temperature**

OrACELL® Decellularized Dermis Allograft

Size	Thickness	Product #	Fee	7 or More
1.0cm x 1.0cm	0.76 - 1.25mm	#DERMIS-1X1	95.00 EACH	84.70 EACH
1.5cm x 2.0cm	0.76 - 1.25mm	#DERMIS-1.5X2	183.00 EACH	164.50 EACH
1.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-1X4	209.00 EACH	187.40 EACH
2.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-2X4	230.00 EACH	206.80 EACH
1.0cm x 1.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X1	95.00 EACH	84.70 EACH
1.5cm x 2.0cm	1.26 - 1.75mm	#DERMIS-THICK-1.5X2	183.00 EACH	164.50 EACH
1.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X4	209.00 EACH	187.40 EACH
2.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-2X4	230.00 EACH	206.80 EACH





OrACELL® Decellularized Dermis Case Studies

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Matt Ficca



Pre-Op
Tooth #6 & #7



6 Year Follow-Up



Pre-Op
Tooth #10 & #11



6 Year Follow-Up



“OrACELL® provides an excellent alternative to autogenous tissue for tissue graft procedures. The pre-hydrated samples handle comparative to autogenous or even better. Over the years I have looked for an allograft with consistent quality and thickness, ease of handling, and predictable healing. I want to offer my patients the highest quality allograft on the market. I am confident I’ve got it covered with OrACELL®.”

—**Dr. Matt Ficca**
Periodontist, Charlotte, NC

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Stephen C. Wallace



Pre-Op



OrACELL® Dermis Placed



Suture Closing




3 Month Post-Op




“Based on the results recorded at 3 months, there was significant gain in root coverage, gain in clinical attachment and increase in width of keratinized gingiva for all recession sites treated with OrACELL® Decellularized Dermis Allograft. I have been using OrACELL® for decades with predictable clinical results.”

—**Dr. Stephen C. Wallace**
Periodontist, Wilmington, NC



Resorbable Bovine Type 1 Collagen Plug

- Use To Stabilize Blood Clots & Preserve Graft Material In Socket Preservation
- Controls Bleeding In Extraction & Biopsy Sites
- Resorbs in 10-14 Days



Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#ORAPLUG-PKG-10	148.00 BOX OF 10	133.50 Box of 10	123.60 Box of 10




Resorbable Collagen Tape

- Use To Repair Sinus Membrane Perforations & Stabilize Blood Clots In Minor Oral Wounds
- Resorbs in 10-14 Days



Resorbs In 10 To 14 Days


Size	Product #	Price	7 or More	14 or More
25mm x 75mm	#ORATAPE-PKG-10	255.00 BOX OF 10	231.70 Box of 10	217.50 Box of 10



80% Calcium Phosphate Granules Combined With 20% Bovine Collagen

- Mineral Mimics Natural Bone In Resorption & Remodeling Profile
- Simple, Cost Effective Socket Preservation
- Plug Takes The Shape Of Defect When Hydrated To Fill Socket
- Offered In Two Sizes - Five Plugs Per Box

Composite Plugs Shown Actual Size - Package Contains 5 Composite Plugs



Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#FULFILL-10X20-PKG-5	275.00 BOX OF 5	258.00 Box of 5	244.30 Box of 5
6mm x 25mm	#FULFILL-6X25-PKG-5	275.00 BOX OF 5	258.00 Box of 5	244.30 Box of 5




Resorbable Collagen Wound Dressings

- For Palatal Donor Sites & To Stabilize Blood Clots In Minor Oral Wounds
- 3X Thicker Than OraTAPE®
- Controls Bleeding & Stabilizes Blood Clots
- Resorbs in 10-14 Days



Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
20mm x 40mm	#ORACOTE-PKG-10	201.00 BOX OF 10	182.70 Box of 10	165.30 Box of 10



Demineralized Cortical Allograft Putty With Cortical Cancellous Chips

- 1-2.8mm Chips Provide An Effective Scaffold For Bone Cell Migration
- Handles Like PlayDoh™ When Hydrated With Blood Or Saline
- Resorbable Hydrogel Carrier Resists Irrigation To Keep Graft In Place




AlloSculpt-3D Case Study


Case & Photos: Dr. Nick Shumaker




Buccal Wall Defect



Bone Chips Provide Scaffold



New Host Bone Formation



20 Week X-Ray Implant Placed



-Dr. Nick Shumaker
Periodontist, Fort Collins, CO

SALVIN AlloSculpt-3D Demineralized Cortical Allograft Putty With Cortical & Cancellous Chips

Size	Product #	Fee	7 or More
0.5 cc	#ALLO-3D-.5	237.00	221.90
1.0 cc	#ALLO-3D-1.0	320.00	295.80
2.0 cc	#ALLO-3D-2.0	470.00	416.60

No Refrigeration Required - Store At Room Temperature
Product Shelf Life - 2.5 Years

Donor Tissue From LifeNet Health® Tissue Bank
Each Lot Is Tested In Vivo For Osteoinductive Potential




PlayDoh™ is a trademark of Hasbro



Scan For

AlloSculpt-3D®
Case Study





A Next Generation Demineralized Bone Matrix (DBM)


- Exceptional Handling – Moldable Putty Conforms To The Graft Site
- Each Lot Tested To Ensure Osteoinductive Potential
- Thickens At Body Temperature, For Better Graft Containment
- Packaged Sterile In A Syringe

Evaluation of Osteoinductive Potential

An In Vitro Study Was Conducted To Examine Accessibility To Bioactive Proteins In ABM And DBM¹. The Content Of BMP-2 Was Measured Over Time Using An Enzyme Linked Immunosorbent Assay (ELISA).

Study Results

- Soluble BMP-2 Was Detected In ABM At An Earlier Time Compared To Particulate DBM
- ABM's Open Pore Structure Provides Earlier Accessibility To Bone Proteins



Accell Connexus




Size	Product #	Fee
0.5 cc	#ACCELL-CONN-0.5CC	332.00
1.0 cc	#ACCELL-CONN-1.0CC	507.00
2.5 cc	#ACCELL-CONN-2.5CC	895.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

1. Khaliq S, Lollis R, Bell D, Oliver R, Walsh WR, and Ingram R, Evaluation of a Next Generation DBM Putty in a Posterolateral Spinal Fusion Model, (2009) Integra LifeSciences Corporation.





Demineralized Cortical Bone With Mineralized Cancellous Chips Paste Or Putty

- For GBR, Ridge Augmentation, Sinus Lifts & Bony Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaBlast Paste



Size	Product #	Fee
0.5 cc	#ORTHO-PASTE-0.5CC	146.00
1.0 cc	#ORTHO-PASTE-1.0CC	203.00
3.0 cc	#ORTHO-PASTE-3.0CC	522.00

DynaBlast Putty



Size	Product #	Fee
1.0 cc	#ORTHO-PUTTY-1CC	249.00
2.5 cc	#ORTHO-PUTTY-2.5CC	498.00
5.0 cc	#ORTHO-PUTTY-5.0CC	881.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

Demineralized Cortical Bone Paste Or Putty

- For Osseous Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaGraft-D Gel



Size	Product #	Fee
0.5 cc	#DYNA-GEL-0.5CC	146.00
1.0 cc	#DYNA-GEL-1.0CC	203.00
5.0 cc	#DYNA-GEL-5.0CC	828.00


DynaGraft-D Putty



Size	Product #	Fee
1.0 cc	#DYNA-PUTTY-1.0CC	249.00
2.5 cc	#DYNA-PUTTY-2.5CC	498.00
5.0 cc	#DYNA-PUTTY-5.0CC	881.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months





Tricalcium Phosphate + Hydroxyapatite

- Promotes Osteoblast Proliferation For New Host Bone Formation
- Alternative To Allografts Or Xenografts For GBR Procedures
- Slow Resorbing, Excellent For Space Maintenance



SALVIN Orastruct 80/20 TCP/HA Synthetic Graft Material

Size	Particulate Size	Product #	Price	7 or More	14 or More
0.5cc	0.5 - 1mm	#ORASTRUCT-0.5CC	86.00	74.70	66.50
1.0cc	0.5 - 1mm	#ORASTRUCT-1.0CC	124.00	108.80	98.40
2.5cc	1 - 2mm	#ORASTRUCT-2.5CC	210.00	191.50	177.60
0.5cc In Syringe - Pkg of 5	0.5 - 1mm	#ORASTRUCT-SYR-0.5CC	309.00	281.70	263.50



Absorbable Collagen Hemostatic Sponges

- Controls Bleeding In 2 To 5 Minutes
- Coherent Sponge Structure Allows Easily Controlled Application
- 0.5" x 1" (1.27cm x 2.54cm)



Package Contains 18 Sponges

Compare To GelFoam®

GelFoam® is a registered trademark of Pharmacia / Upjohn Company

SALVIN OraFOAM Absorbable Collagen Hemostatic Sponges

Size	Product #	Price	2 or More Boxes
0.5" x 1" (1.27cm x 2.54cm)	#ORAFOAM-PKG-18	311.00 BOX OF 18	296.30 Box of 18

OraFOAM™ Hemostatic Sponge Is Indicated In Surgical Procedures (Other Than Ophthalmological & Urological Surgery) As An Adjunct To Hemostasis When Control Of Bleeding By Ligature Or Conventional Procedures Is Ineffective Or Impractical. Please See Product IFU For Precautions & Adverse Reactions. IFU Available Online At www.salvin.com/Downloads/Guides-Studies.html



Non-Resorbable PTFE Sutures

- Excellent Choice For Bone Grafting & Implant Procedures
- Smooth & Flexible – Not Fragile Or Brittle – So Knots Stay Tied
- Proprietary Monofilament Is Biologically Inert & Non-Cytotoxic
- Longer Length – All PTFE Sutures Now 24" Long



Black Needle For Higher Visibility & Less Glare


Super Sharp Needle Requires Less Force – Easier Suturing With Less Trauma

"Racetrack" Rather Than Zig-Zag Package




Little-To-No Package Memory For Excellent Handling And Secure Knotting


Ultra-Soft, Ultra-Smooth Suture Designed With Low Surface Porosity



SALVIN CytoSurg



Competitor G (2000x Magnification)




Competitor O

	Salvin Item#	Size	Needle	Length	Per Box	10 Or More
	#PTFE-20FS2	2-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	112.60 Box 12	100.90 Box 12
MOST POPULAR	#PTFE-30FS2	3-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	112.60 Box 12	100.90 Box 12
	#PTFE-30FS3	3-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	112.60 Box 12	100.90 Box 12
	#PTFE-40TF	4-0	13mm 1/2 Circle - <u>Taper Point</u>	24" (60cm)	112.60 Box 12	100.90 Box 12
MOST POPULAR	#PTFE-40FS3	4-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	112.60 Box 12	100.90 Box 12
	#PTFE-40PS2	4-0	19mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	133.20 Box 12	119.50 Box 12
	#PTFE-40P3	4-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	133.20 Box 12	119.50 Box 12
	#PTFE-50P3	5-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	133.20 Box 12	119.50 Box 12
	#PTFE-60P3	6-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	133.20 Box 12	119.50 Box 12

**For Our Full Line of Sutures,
Check Out The
2024 Salvin Instrument Catalog**










SALVIN® SurgiPoint® PGA Polyglycolic Acid Synthetic Absorbable-Braided



- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl™ But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 21-28 Days





	Salvin Item#	Size	Needle	Length	Color	Vicryl™ Comparison	Per Box	10 Or More
▽  FS-2 (C6)	#LV-423B	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J423G	85.50 Box 12	80.70 Box 12
	#LV-422B	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J422G	85.50 Box 12	80.70 Box 12
	#LV-MC397NC	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Violet		125.60 Box 12	118.10 Box 12
	#LV-421B	5-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J421G	85.50 Box 12	80.70 Box 12
▽  FS-3 (C22)	#LV-500B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	White	J500G	81.90 Box 12	78.10 Box 12
▽  P-3 (C3)	#LV-493B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J493G	148.50 Box 12	140.00 Box 12
	#LV-492B	6-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J492G	148.50 Box 12	140.00 Box 12
▽  PS-2 (PC31)	#LV-497B	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J497G	148.50 Box 12	140.00 Box 12
	#LV-496B	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J496G	148.50 Box 12	140.00 Box 12
	#LV-495B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J495G	148.50 Box 12	140.00 Box 12
○  SH-1 (T14)	#LV-311B	3-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J311G	76.20 Box 12	70.90 Box 12
	#LV-310B	4-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J310G	76.20 Box 12	70.90 Box 12

SALVIN® SurgiPoint® PGA Polyglycolic Acid Synthetic Absorbable-Braided Fast Absorbing



- Good Alternative To Plain Gut Sutures
- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl Rapide™ But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 7-14 Days

	Salvin Item#	Size	Needle	Length	Color	Vicryl™ Comparison	Per Box	10 Or More
▽  FS-2 (C6)	#LVFA-493N	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VR497	100.10 Box 12	94.10 Box 12
	#LVFA-423N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed		100.10 Box 12	94.10 Box 12
	#LVFA-492N	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VC422H	100.10 Box 12	94.10 Box 12
	#LVFA-422N	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR426	100.10 Box 12	94.10 Box 12
▽  FS-1 (C7)	#LVFA-442N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR935	100.10 Box 12	94.10 Box 12
▽  PS-2 (PC31)	#LVFA-497N	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		100.10 Box 12	94.10 Box 12
	#LVFA-496N	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		100.10 Box 12	94.10 Box 12
▽  P-1 (C1)	#LVFA-490C	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		129.60 Box 12	119.40 Box 12

More Sutures On Next 2 Pages →

SALVIN® SurgiPoint® Silk **Non-Resorbable**



• Non-Resorbable, Braided, Coated Silk – Excellent Handling & Tying Characteristics

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS (C8)	#L-770B	2-0	3/8 Circle - Reverse Cutting	18" (45cm)	685-G	56.60 Box 12	52.70 Box 12
FS-2 (C6) MOST POPULAR	#L-784B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		56.60 Box 12	52.70 Box 12
	#L-781B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	683G	56.60 Box 12	52.70 Box 12
	#L-754B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	682G	56.60 Box 12	52.70 Box 12
FS-1 (C7)	#L-786B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	684G	56.60 Box 12	52.70 Box 12
C-2 (C17)	#L-787B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	736G	56.60 Box 12	52.70 Box 12
X-1 (C31)	#L-785B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	632G	56.60 Box 12	52.70 Box 12

SALVIN® SurgiPoint® Polypropylene **Non-Resorbable**



• Non-Resorbable Blue Monofilament – Inert, Easy To Handle, Passes Easily Through Tissue

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
P-3 (C3) MOST POPULAR	#L-1085B	5-0	3/8 Circle - Premium Reverse Cutting	18" (45cm)	8698G	164.90 Box 12	153.60 Box 12
	#L-1082B	6-0	3/8 Circle - Premium Reverse Cutting	18" (45cm)	8695G	164.90 Box 12	153.60 Box 12
P-1 (C1)	#L-1079B	6-0	3/8 Circle - Premium Reverse Cutting	18" (45cm)	8697G	164.90 Box 12	153.60 Box 12
	#L-2468N	7-0	3/8 Circle - Premium Reverse Cutting	18" (45cm)	8696G	164.90 Box 12	153.60 Box 12
C-2 (C17)	#L-1011B	6-0	3/8 Circle - Reverse Cutting	10" (25cm)	8889H	164.90 Box 12	153.60 Box 12
FS-2 (C6)	#L-8661B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8661G	84.90 Box 12	80.90 Box 12
FS-3 (C22)	#L-8680B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8680G	84.90 Box 12	80.90 Box 12

SALVIN® SurgiPoint® Nylon Monofilament **Non-Resorbable**



• Non-Resorbable Black Monofilament – Passes Easily Through Tissue – Similar To Resolon™


	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-922B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	662G	65.90 Box 12	61.60 Box 12
	#L-918B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	661G	65.90 Box 12	61.60 Box 12
P-3 (C3)	#L-698BKC	5-0	3/8 Circle - Black Premium Reverse	18" (45cm)	698G	117.20 Box 12	111.70 Box 12
	#L-1698B	6-0	3/8 Circle - Premium Reverse Cutting	18" (45cm)	1698G	117.20 Box 12	111.70 Box 12

SALVIN® SurgiPoint® Polyviolene-Green Braided Polyester **Non-Resorbable**



• Green Braided Uncoated Polyester – Smooth, Easy To Handle With Good Knot Security






	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-1180B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	R683G	56.60 Box 12	52.70 Box 12



SALVIN® SurgiPoint® Chromic Gut *Resorbs In 7-14 Days*







- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6) <small>MOST POPULAR</small>	#L-560B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	636G	56.70 Box 12	52.80 Box 12
	#L-559B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	635G	56.70 Box 12	52.80 Box 12
	#L-525B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		56.70 Box 12	52.80 Box 12
	#L-558B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		56.70 Box 12	52.80 Box 12
	#L-557B	5-0	3/8 Circle - Reverse Cutting	27" (68cm)		56.70 Box 12	52.80 Box 12
 P-3 (C3) <small>MOST POPULAR</small>	#L-595B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1654G	107.60 Box 12	101.60 Box 12
	#L-1248B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	687G	107.60 Box 12	101.60 Box 12
 P-1 (C1)	#1242B-CP	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		180.30 Box 12	169.20 Box 12
 SH-1 (T4)	#L-589B	4-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		69.40 Box 12	64.50 Box 12
 RB-1 (T28)	#L-590B	5-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		56.70 Box 12	53.90 Box 12

SALVIN® SurgiPoint® Gut-Plain *Resorbs In 5-10 Days*







- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6) <small>MOST POPULAR</small>	#L-554B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	56.70 Box 12	52.80 Box 12
	#L-551B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	56.70 Box 12	52.80 Box 12
	#L-591B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		56.70 Box 12	52.80 Box 12
	#L-592B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		56.70 Box 12	52.80 Box 12
 P-3 (C3)	#L-596B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	686G	107.60 Box 12	101.60 Box 12
 X-1 (C31)	#L-594B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	612G	56.70 Box 12	52.80 Box 12
 PC-1 (CP1)	#L-1415B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1916G	182.30 Box 12	171.10 Box 12

SALVIN® SurgiPoint® Monoderm™ *Resorbs In 21-28 Days*



- Good Alternative To Chromic Gut Sutures
- Synthetic Resorbable Monofilament Glycolide & E-Caprolactone – Similar to Monocryl™

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6)	#L-388	3-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		90.60 Box 12	82.50 Box 12
	#L-392C	4-0	3/8 Circle - Black Premium Reverse	18" (45cm)		90.60 Box 12	82.50 Box 12
 FS-3 (C22)	#L-395	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		90.60 Box 12	82.50 Box 12
 PS-2 (DGL19) <small>MOST POPULAR</small>	L-YD513N	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	Y513G	145.70 Box 12	137.70 Box 12
 P-3 (C3) <small>NEW!</small>	#L-464B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		152.60 Box 12	144.40 Box 12
	#L-463B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		152.60 Box 12	144.40 Box 12
	#L-463C	5-0	3/8 Circle - Black Premium Reverse	18" (45cm)	Y463G	152.60 Box 12	144.40 Box 12
	#L-492B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		152.60 Box 12	144.40 Box 12
	#L-Q499B	6-0	3/8 Circle - Violet Premium Reverse	18" (45cm)		152.60 Box 12	144.40 Box 12

Suture Names Are Trademarks Of Their Respective Companies.



Toll Free 800.535.6566
www.salvin.com

Salvin Dental Specialties, Inc.

3450 Latrobe Drive
Charlotte, NC 28211 U.S.A.

CHANGE SERVICE REQUESTED

PRSRRT STD
U.S. POSTAGE
PAID
PPCO

**Whether You Are New To Implant Dentistry Or Performing
Advanced Procedures, We've Got The Solutions You Need To
Make Your Life Easier And Your Practice More Productive.**

SALVIN[®]
Dental Specialties

SALVIN[®]
Regenerative

***Everything For Your Implant Practice
But The Implants[®]***